


Npy2r Antibody (R34349)

Catalog No.	Formulation	Size
R34349-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4887
Applications	Western Blot : 0.1-0.3ug/ml Immunoprecipitation : suitable ELISA (peptide) LOD : 1:32000
Limitations	This Npy2r antibody is available for research use only.



Western blot testing of rat brain lysate with Npy2r antibody at 0.1ug/ml. Predicted molecular weight: ~43 kDa.

Description

Additional name(s) for this target protein: Neuropeptide Y receptor Y2

Application Notes

Optimal dilution of the Npy2r antibody should be determined by the researcher.

Immunogen

Amino acids KAKKNLEVKKNN were used as the immunogen for this Npy2r antibody.

Storage

Aliquot and store the Npy2r antibody at -20oC.