


NPY5R Antibody Biotin Conjugate (R34109BTN)

Catalog No.	Formulation	Size
R34109BTN	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4889
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:32000
Limitations	This NPY5R antibody is available for research use only.



Description

Additional name(s) for this target protein: Neuropeptide Y receptor Y5

Application Notes

Optimal dilution of the NPY5R antibody should be determined by the researcher.

Immunogen

Amino acids GLSNKENRLEEN were used as the immunogen for this NPY5R antibody.

Storage

Aliquot and store the NPY5R antibody at -20oC.