

# NPY5R Antibody (R34109)


Catalog No.	Formulation	Size
R34109-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4889
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : 4-6ug/ml Immunoprecipitation : suitable ELISA (peptide) LOD : 1:32000
Limitations	This NPY5R antibody is available for research use only.

KUa
200 -
116 –
97 –
66 –
44 –
31 —
22 –
14 —
6

Western blot testing of human brain lysate with NPY5R antibody at 0.3ug/ml. Predicted/observed molecular weight: ~51kDa.



IHC staining of FFPE human intestine with NPY5R antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

# **Description**

Additional name(s) for this target protein: Neuropeptide Y receptor Y5

# **Application Notes**

Optimal dilution of the NPY5R antibody should be determined by the researcher.

### **Immunogen**

Amino acids GLSNKENRLEEN were used as the immunogen for this NPY5R antibody.

#### **Storage**

Aliquot and store the NPY5R antibody at -20oC.