


NR0B1 Antibody / DAX-1 (R33235)

Catalog No.	Formulation	Size
R33235-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	190
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This NR0B1 antibody is available for research use only.



Western blot testing of human testis lysate with NR0B1 antibody at 0.5ug/ml. Predicted molecular weight: 44-52 kDa (two isoforms).

Description

Additional name(s) for this target protein: Nuclear receptor subfamily 0, group B, member 1; DAX1; Dosage-sensitive sex reversal, adrenal hypoplasia critical region, on chromosome X, gene 1

Application Notes

Optimal dilution of the NR0B1 antibody should be determined by the researcher.

Immunogen

Amino acids CGEDHPQQGSTLY were used as the immunogen for this NR0B1 antibody.

Storage

Aliquot and store the NR0B1 antibody at -20oC.