

NRAS Antibody (R35682)


Catalog No.	Formulation	Size
R35682-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4893
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 2-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This NRAS antibody is available for research use only.



Western blot testing of human A431 lysate with NRAS antibody at 0.01ug/ml. Predicted molecular weight ~21 kDa.



Western blot testing of rat intestine lysate with NRAS antibody at 0.1ug/ml. Predicted molecular weight ~21 kDa.



IHC staining of FFPE human small intestine with NRAS antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Neuroblastoma RAS

Application Notes

Optimal dilution of the NRAS antibody should be determined by the researcher.

Immunogen

Amino acids RMKKLNSSDDGTQ were used as the immunogen for this NRAS antibody.

Storage

Aliquot and store the NRAS antibody at -20oC.