


## NRG3 Antibody (R35689)

Catalog No.	Formulation	Size
R35689-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10718
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This NRG3 antibody is available for research use only.



Western blot testing of human cerebellum lysate with NRG3 antibody at 0.3ug/ml. The expected ~75kDa band and the additional ~37kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: Neuregulin 3

## Application Notes

Optimal dilution of the NRG3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids RNEIQRDSALTK were used as the immunogen for this NRG3 antibody.

## **Storage**

Aliquot and store the NRG3 antibody at -20oC.