

Nsg2 Antibody / Neuron Specific Gene 2 (R36446)

Catalog No.	Formulation	Size
R36446-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P47759
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) : 1:128000
Limitations	This Nsg2 antibody is available for research use only.



Western blot testing of mouse fetal brain lysate with Nsg2 antibody at 0.1ug/ml.
Predicted molecular weight ~19 kDa.

Description

Additional name(s) for this target protein: Neuron specific gene family member 2

Application Notes

The stated application concentrations are suggested starting points. Titration of the Nsg2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Amino acids HEPKPPKTQGH were used as the immunogen for this Nsg2 antibody.

Storage

Aliquot and store the Nsg2 antibody at -20oC.