

## Nsg2 Antibody / Neuron Specific Gene 2 (R36446)

Catalog No.	Formulation	Size
R36446-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Mouse
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	P47759
<b>Applications</b>	Western Blot : 0.1-0.3ug/ml ELISA (peptide) : 1:128000
<b>Limitations</b>	This Nsg2 antibody is available for research use only.



Western blot testing of mouse fetal brain lysate with Nsg2 antibody at 0.1ug/ml.  
Predicted molecular weight ~19 kDa.

### Description

Additional name(s) for this target protein: Neuron specific gene family member 2

### Application Notes

The stated application concentrations are suggested starting points. Titration of the Nsg2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## **Immunogen**

Amino acids HEPKPPKTQGH were used as the immunogen for this Nsg2 antibody.

## **Storage**

Aliquot and store the Nsg2 antibody at -20oC.