


## NT5C3 Antibody (F53927)

| Catalog No.   | Formulation                                | Size    |
|---------------|--------------------------------------------|---------|
| F53927-0.2ML  | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.2 ml  |
| F53927-0.05ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.05 ml |

[Bulk quote request](#)

|                             |                                                         |
|-----------------------------|---------------------------------------------------------|
| <b>Availability</b>         | 1-3 business days                                       |
| <b>Species Reactivity</b>   | Human                                                   |
| <b>Predicted Reactivity</b> | Mouse                                                   |
| <b>Format</b>               | Antigen affinity purified                               |
| <b>Host</b>                 | Rabbit                                                  |
| <b>Clonality</b>            | Polyclonal (rabbit origin)                              |
| <b>Isotype</b>              | Rabbit Ig                                               |
| <b>Purity</b>               | Antigen affinity                                        |
| <b>UniProt</b>              | Q9H0P0                                                  |
| <b>Applications</b>         | Western Blot : 1:1000-2000                              |
| <b>Limitations</b>          | This NT5C3 antibody is available for research use only. |



Western blot testing of human 1) A431, 2) HeLa, 3) HepG2, 4) Jurkat, 5) U-2OS and 6) SH-SY5Y cell lysate with NT5C3 antibody at 1:2000. Predicted molecular weight: 38 kDa.

## Description

Cytosolic 5'-nucleotidase 3A can act both as nucleotidase and as phosphotransferase.

## Application Notes

Titration of the NT5C3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 282-313 from the human protein was used as the immunogen for the NT5C3 antibody.

## Storage

Aliquot the NT5C3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.