

NTRK2 Antibody / Neurotrophic receptor tyrosine kinase 2 [clone NTRK2/7925] (V5960)

| Catalog No. | Formulation | Size |
|----------------|--|--------|
| V5960-100UG | 0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide | 100 ug |
| V5960-20UG | 0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide | 20 ug |
| V5960SAF-100UG | 1 mg/ml in 1X PBS; BSA free, sodium azide free | 100 ug |

[Bulk quote request](#)

| | |
|---------------------------|---|
| Species Reactivity | Human |
| Format | Purified |
| Host | Mouse |
| Clonality | Monoclonal (mouse origin) |
| Isotype | Mouse IgG2, kappa |
| Clone Name | NTRK2/7925 |
| UniProt | Q16620 |
| Localization | Cell surface. |
| Applications | Western Blot : 2-4ug/ml |
| Limitations | This NTRK2/Neurotrophic receptor tyrosine kinase 2 antibody is available for research use only. |

Application Notes

Optimal dilution of the NTRK2/Neurotrophic receptor tyrosine kinase 2 antibody should be determined by the researcher.

Immunogen

A recombinant fragment (around amino acids 250-450) of human NTRK2 protein (exact sequence is proprietary) was used as the immunogen for the NTRK2/Neurotrophic receptor tyrosine kinase 2 antibody.

Storage

NTRK2/Neurotrophic receptor tyrosine kinase 2 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

