


NTRK2 Antibody (R34672)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34672-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Mouse, Rat, Dog |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 4915 |
| Applications | Western Blot : 0.05-0.2ug/ml ELISA (peptide) LOD : 1:64000 |
| Limitations | This NTRK2 antibody is available for research use only. |



Western blot testing of human hippocampus [A], cerebral cortex [B] and cerebellum [C] lysates with NTRK2 antibody at 0.05ug/ml. Expected molecular weight: 90~140 kDa, depending on glycosylation level.

Description

Additional name(s) for this target protein: Neurotrophic tyrosine kinase, receptor, type 2

Application Notes

Optimal dilution of the NTRK2 antibody should be determined by the researcher.

1. This NTRK2 antibody will detect isoforms a/b/c/d/e.

Immunogen

Amino acids KTLQEAKSSPDTQ were used as the immunogen for this NTRK2 antibody.

Storage

Aliquot and store the NTRK2 antibody at -20oC.