


NUDT19 Antibody / Nudix motif 19 [clone 29N64] (RQ8688)

Catalog No.	Formulation	Size
RQ8688	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	29N64
Purity	Affinity chromatography
UniProt	A8MXV4
Applications	Western Blot : 1:500-1:2000
Limitations	This NUDT19 antibody is available for research use only.



Western blot testing of human 293T cell lysate with NUDT19 antibody. Predicted molecular weight ~42 kDa.

Description

NUDT19 is an acyl-CoA diphosphatase that mediates the hydrolysis of a wide range of CoA esters, including choloyl-CoA and branched-chain fatty-acyl-CoA esters.

Application Notes

Optimal dilution of the NUDT19 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human NUDT19 protein was used as the immunogen for the NUDT19 Antibody.

Storage

Store the NUDT19 antibody at -20oC.