


NUMBL Antibody (R36056)

Catalog No.	Formulation	Size
R36056-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9253
Applications	Western Blot : 0.5-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This NUMBL antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human NUMBL-MYC probed with NUMBL antibody [1ug/ml] in Lane A and probed with anti-MYC [1/1000] in Lane C. Mock-transfected HEK293 probed with NUMBL antibody [1ug/ml] in Lane B. Observed molecular weight ~85kDa.

Description

Additional name(s) for this target protein: Numb-like

Application Notes

Optimal dilution of the NUMBL antibody should be determined by the researcher.

Immunogen

Amino acids SRSAAASGGPRRPE were used as the immunogen for this NUMBL antibody. **Storage** Aliquot and store the NUMBL antibody at -20oC.