


## NUMBL Antibody (R36056)

Catalog No.	Formulation	Size
R36056-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9253
Applications	Western Blot : 0.5-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This NUMBL antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human NUMBL-MYC probed with NUMBL antibody [1ug/ml] in Lane A and probed with anti-MYC [1/1000] in Lane C. Mock-transfected HEK293 probed with NUMBL antibody [1ug/ml] in Lane B. Observed molecular weight ~85kDa.

## Description

Additional name(s) for this target protein: Numb-like

## Application Notes

Optimal dilution of the NUMBL antibody should be determined by the researcher.

## Immunogen

Amino acids SRSAAASGGPRRPE were used as the immunogen for this NUMBL antibody.

## Storage

Aliquot and store the NUMBL antibody at -20oC.