


NUP107 Antibody (R35062)


Catalog No.	Formulation	Size
R35062-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57122
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This NUP107 antibody is available for research use only.



Western blot testing of peripheral blood lymphocyte lysate with NUP107 antibody at 1ug/ml. Expected molecular weight: ~107kDa.



Western blot testing of nuclear lysate of NIH3T3 with NUP107 antibody at 2ug/ml. Expected molecular weight: ~107kDa.

Description

Additional name(s) for this target protein: Nucleoporin 107kDa, NUP84

Application Notes

Optimal dilution of the NUP107 antibody should be determined by the researcher.

Immunogen

Amino acids PLPAEDDNAIREH were used as the immunogen for this NUP107 antibody.

Storage

Aliquot and store the NUP107 antibody at -20oC.