


Nup50 Antibody (R33863)


Catalog No.	Formulation	Size
R33863-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10762
Applications	Western Blot : 0.2-1ug/ml Immunohistochemistry (FFPE) : 2.5-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Nup50 antibody is available for research use only.



Western blot testing of Jurkat lysate with Nup50 antibody at 0.5ug/ml. Predicted molecular weight: observed at ~50 or ~60kDa.



IHC staining of FFPE human uterus with Nup50 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Nucleoporin 50kDa; Npap60

Application Notes

Optimal dilution of the Nup50 antibody should be determined by the researcher.

Immunogen

Amino acids ELHKILLEKKDA were used as the immunogen for this Nup50 antibody.

Storage

Aliquot and store the Nup50 antibody at -20oC.