


OAS1 Antibody (R36183)

Catalog No.	Formulation	Size
R36183-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4938
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This OAS1 antibody is available for research use only.



Western blot testing of HepG2 lysate with OAS1 antibody at 1ug/ml. Predicted molecular weight: ~46kDa.

Description

Additional name(s) for this target protein: 2'-5'-oligoadenylate synthetase 1

Application Notes

Optimal dilution of the OAS1 antibody should be determined by the researcher.

1. This OAS1 antibody is specific for isoform 1.

Immunogen

Amino acids THEYPHFSHRPST were used as the immunogen for this OAS1 antibody.

Storage

Aliquot and store the OAS1 antibody at -20oC.