


OAS2 Antibody (R35659)

Catalog No.	Formulation	Size
R35659-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4939
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This OAS2 antibody is available for research use only.



Western blot testing of Jurkat lysate with OAS2 antibody at 1ug/ml. Image shows bands at ~85kDa and ~19kDa.

Description

Additional name(s) for this target protein: 2'-5'-oligoadenylate synthase 2

Application Notes

Optimal dilution of the OAS2 antibody should be determined by the researcher.

Immunogen

Amino acids RKTVLRGNSDGT were used as the immunogen for this OAS2 antibody. This sequence is common to isoforms 1 (82 kDa, can be observed at 71kDa) and 3 (~20kDa).

Storage

Aliquot and store the OAS2 antibody at -20oC.