


OAS3 Antibody (R34249)

Catalog No.	Formulation	Size
R34249-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4940
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This OAS3 antibody is available for research use only.



Western blot testing of A431 lysate with OAS3 antibody at 2ug/ml. Predicted molecular weight: ~121kDa, observed here at ~110kDa.

Description

Additional name(s) for this target protein: 2'-5'-oligoadenylate synthetase 3

Application Notes

Optimal dilution of the OAS3 antibody should be determined by the researcher.

Immunogen

Amino acids HPIQLDPNQKTPENS were used as the immunogen for this OAS3 antibody.

Storage

Aliquot and store the OAS3 antibody at -20oC.