


## OAS3 Antibody (R34249)

Catalog No.	Formulation	Size
R34249-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4940
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This OAS3 antibody is available for research use only.



Western blot testing of A431 lysate with OAS3 antibody at 2ug/ml. Predicted molecular weight: ~121kDa, observed here at ~110kDa.

## Description

Additional name(s) for this target protein: 2'-5'-oligoadenylate synthetase 3

## Application Notes

Optimal dilution of the OAS3 antibody should be determined by the researcher.

## Immunogen

Amino acids HPIQLDPNQKTPENS were used as the immunogen for this OAS3 antibody.

## Storage

Aliquot and store the OAS3 antibody at -20oC.