

## OAZ1 Antibody / Ornithine decarboxylase antizyme 1 (R36432)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R36432-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |   |
|--------------------|---|
| Availability       | 1-3 business days   |
| Species Reactivity | Human   |
| Format             | Antigen affinity purified                                   |
| Clonality          | Polyclonal (goat origin)                                    |
| Isotype            | Goat Ig   |
| Purity             | Antigen affinity  |
| Gene ID            | 4946  |
| Applications       | Western Blot : 0.3-1ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations        | This OAZ1 antibody is available for research use only.      |



Western blot testing of human kidney lysate with OAZ1 antibody at 0.3ug/ml. Expected molecular weight: ~25 kDa.

### Description

Additional name(s) for this target protein: OAZ, AZI, ODC-Az

### Application Notes

Optimal dilution of the OAZ1 antibody should be determined by the researcher.

### Immunogen

Amino acids SRLTDAKRINWRT were used as the immunogen for this OAZ1 antibody.

## Storage

Aliquot and store the OAZ1 antibody at -20oC.