

OAZ1 Antibody / Ornithine decarboxylase antizyme 1 (R36432)

Catalog No.	Formulation	Size
R36432-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4946
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This OAZ1 antibody is available for research use only.



Western blot testing of human kidney lysate with OAZ1 antibody at 0.3ug/ml. Expected molecular weight: ~25 kDa.

Description

Additional name(s) for this target protein: OAZ, AZI, ODC-Az

Application Notes

Optimal dilution of the OAZ1 antibody should be determined by the researcher.

Immunogen

Amino acids SRLTDAKRINWRT were used as the immunogen for this OAZ1 antibody.

Storage

Aliquot and store the OAZ1 antibody at -20oC.