


OCT-2 Antibody / POU2F2 [clone AEBG-16] (RQ5266)

Catalog No.	Formulation	Size
RQ5266	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AEBG-16
Purity	Affinity purified
UniProt	P09086
Applications	Western Blot : 1:500-1:2000
Limitations	This OCT-2 antibody is available for research use only.



Western blot testing with OCT-2 antibody. Expected molecular weight: isoforms from 43~62 kDa.

Description

Transcription factor that specifically binds to the octamer motif (5'-ATTTGCAT-3'). Regulates transcription in a number of tissues in addition to activating immunoglobulin gene expression.

Application Notes

Optimal dilution of the OCT-2 antibody should be determined by the researcher.

Immunogen

Full length human OCT-2 / POU2F2 protein was used as the immunogen for the OCT-2 antibody.

Storage

Store the OCT-2 antibody at -20oC.