


OCT6 Antibody (R34171)

Catalog No.	Formulation	Size
R34171-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5453
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This OCT6 antibody is available for research use only.



Western blot testing of A549 lysate with OCT6 antibody at 1ug/ml. Predicted/observed molecular weight: ~45kDa.

Description

Additional name(s) for this target protein: Octamer-binding transcription factor 6; POU3F1

Application Notes

Optimal dilution of the OCT6 antibody should be determined by the researcher.

Immunogen

Amino acids HEDGHEAQLEPSP were used as the immunogen for this OCT6 antibody.

Storage

Aliquot and store the OCT6 antibody at -20oC.