


ODF2 Antibody (R35590)

Catalog No.	Formulation	Size
R35590-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4957
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ODF2 antibody is available for research use only.



Western blot testing of C57/B6 mouse CTL lysate with ODF2 antibody at 3ug/ml.
Predicted molecular weight: ~95kDa, routinely observed at ~84kDa.

Description

Additional name(s) for this target protein: Outer dense fiber of sperm tails 2, 84 kD; ODF84

Application Notes

Optimal dilution of the ODF2 antibody should be determined by the researcher.

1. This ODF2 antibody is specific for isoform 1.

Immunogen

Amino acids REKHQASQKENKQ were used as the immunogen for this ODF2 antibody.

Storage

Aliquot and store the ODF2 antibody at -20oC.