


OLIG2 Antibody [clone OLIG2/2400] (V3819)

Catalog No.	Formulation	Size
V3819-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V3819-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V3819SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug


[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	OLIG2/2400
Purity	Protein G affinity
UniProt	Q13516
Localization	Nuclear
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT
Limitations	This OLIG2 antibody is available for research use only.




IHC testing of FFPE human cerebellum with OLIG2 antibody (clone OLIG2/2400). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 10-20 min followed by cooling at RT for 20 min.

Human Protein Microarray Specificity Validation



Analysis of HuProt(TM) microarray containing more than 19,000 full-length human proteins using OLIG2 antibody (clone OLIG2/2400). These results demonstrate the foremost specificity of the OLIG2/2400 mAb.



SDS-PAGE analysis of purified, BSA-free OLIG2 antibody (clone OLIG2/2400) as confirmation of integrity and purity.

Description

Olig2, a basic helix-loop-helix transcription factor, is involved in oligodendroglial specification. Olig2 expression has been reported in most glial tumors, such as oligodendrogliomas and astrocytomas. Although more than half of glioblastomas are positive for Olig2, expression is very weak in terms of both percentage of labeled cells and intensity. No Olig2 expression has been found in the non-glial tumors including neuro-epithelial tumors, ependymomas, sub-ependymomas, medulloblastomas, and non-neuroepithelial tumors, such as CNS lymphomas, meningiomas, schwannomas, atypical teratoid / rhabdoid tumor, and haemangioblastomas. Compared to the strong staining seen in glioma samples, a weak expression is observed in non-tumoral brain tissue (gliosis).

Application Notes

Optimal dilution of the OLIG2 antibody should be determined by the researcher.

Immunogen

A portion of amino acids 1-141 from the human protein was used as the immunogen for this OLIG2 antibody.

Storage

After reconstitution, the OLIG2 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.