


## Olig3 Antibody (F53071)

Catalog No.	Formulation	Size
F53071-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53071-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q7RTU3
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This Olig3 antibody is available for research use only.



Western blot analysis of lysate from (1) mouse brain and (2) mouse skeletal muscle tissue lysate using Olig3 antibody at 1:1000.

## Description

May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons.

## Application Notes

Titration of the Olig3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## **Immunogen**

This Olig3 antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 15-48 amino acids from the N-terminal region of human OLIG3.

## **Storage**

Aliquot the Olig3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.