


ORC4L Antibody (R33669)


Catalog No.	Formulation	Size
R33669-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5000
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This ORC4L antibody is available for research use only.



Western blot testing of mouse heart lysate with ORC4L antibody at 0.5ug/ml. Predicted molecular weight: ~50 kDa.



Western blot testing of Jurkat nuclear lysate with ORC4L antibody at 0.5ug/ml. Predicted molecular weight: ~50 kDa.

Description

Additional name(s) for this target protein: Origin recognition complex, subunit 4-like

Application Notes

Optimal dilution of the ORC4L antibody should be determined by the researcher.

Immunogen

Amino acids DVRQWATSSLSWL were used as the immunogen for this ORC4L antibody.

Storage

Aliquot and store the ORC4L antibody at -20oC.