


## OSBP Antibody (R34462)

Catalog No.	Formulation	Size
R34462-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5007
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This OSBP antibody is available for research use only.



Western blot testing of human placenta lysate with OSBP antibody at 0.1ug/ml.  
Predicted molecular weight: ~89kDa.

## Description

Additional name(s) for this target protein: Oxysterol binding protein

## Application Notes

Optimal dilution of the OSBP antibody should be determined by the researcher.

## Immunogen

Amino acids KHQLEETKKEKRTRIP were used as the immunogen for this OSBP antibody.

## Storage

Aliquot and store the OSBP antibody at -20oC.