


OSBP1 Antibody (R33260)

Catalog No.	Formulation	Size
R33260-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5007
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This OSBP1 antibody is available for research use only.



Western blot testing of human liver lysate with OSBP1 antibody at 2ug/ml. Predicted molecular weight: ~90kDa.

Description

Additional name(s) for this target protein: Oxysterol binding protein; OSBP

Application Notes

Optimal dilution of the OSBP1 antibody should be determined by the researcher.

Immunogen

Amino acids CKEKQDWSSCPDIF were used as the immunogen for this OSBP1 antibody.

Storage

Aliquot and store the OSBP1 antibody at -20oC.