


OSBPL10 Antibody (R33369)

Catalog No.	Formulation	Size
R33369-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	114884
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This OSBPL10 antibody is available for research use only.



Western blot testing of human duodenum lysate with OSBPL10 antibody at 1ug/ml.
Predicted molecular weight: ~84kDa.

Description

Additional name(s) for this target protein: Oxysterol binding protein-like 10

Application Notes

Optimal dilution of the OSBPL10 antibody should be determined by the researcher.

Immunogen

Amino acids DAATEQKRHLEEKQ were used as the immunogen for this OSBPL10 antibody.

Storage

Aliquot and store the OSBPL10 antibody at -20oC.