


OSBPL11 Antibody (R34570)

Catalog No.	Formulation	Size
R34570-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	114885
Applications	Western Blot : 1-2ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This OSBPL11 antibody is available for research use only.



Western blot testing of human brain lysate with OSBPL11 antibody at 1ug/ml. Predicted molecular weight: ~86kDa.

Description

Additional name(s) for this target protein: Oxysterol-binding protein-related protein 11

Application Notes

Optimal dilution of the OSBPL11 antibody should be determined by the researcher.

Immunogen

Amino acids KPLWKIIPTTQPAE were used as the immunogen for this OSBPL11 antibody.

Storage

Aliquot and store the OSBPL11 antibody at -20oC.