


OSBPL2 Antibody (R36366)

Catalog No.	Formulation	Size
R36366-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9885
Applications	Western Blot : 1-2ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This OSBPL2 antibody is available for research use only.



Western blot testing of human heart lysate with OSBPL2 antibody at 1.5ug/ml. Predicted molecular weight: ~55kDa.

Description

Additional name(s) for this target protein: Oxysterol binding protein-like 2; ORP2; ORP-2

Application Notes

Optimal dilution of the OSBPL2 antibody should be determined by the researcher.

Immunogen

Amino acids YFERNFSDCPDIY were used as the immunogen for this OSBPL2 antibody.

Storage

Aliquot and store the OSBPL2 antibody at -20oC.