

OSBPL5 Antibody (R34396)


Catalog No.	Formulation	Size
R34396-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	114879
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This OSBPL5 antibody is available for research use only.



Western blot testing of human spleen lysate with OSBPL5 antibody at 2ug/ml. Predicted molecular weight: ~98 kDa, observed here at ~110 kDa.



Western blot testing of human spleen lysate with OSBPL5 antibody at 0.2ug/ml. Predicted molecular weight: ~98 kDa, observed here at ~110 kDa.

Description

Additional name(s) for this target protein: Oxysterol-binding protein-like protein 5, OSBPL5, ORP5, OBPH1

Application Notes

Optimal dilution of the OSBPL5 antibody should be determined by the researcher.

Immunogen

Amino acids KEEAFLRRRFSLC were used as the immunogen for this OSBPL5 antibody.

Storage

Aliquot and store the OSBPL5 antibody at -20oC.