


Osteopontin Antibody (R33573)

Catalog No.	Formulation	Size
R33573-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6696
Localization	Cytoplasmic
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Osteopontin antibody is available for research use only.



Description

Additional name(s) for this target protein: Secreted phosphoprotein 1; SPP1

Application Notes

Optimal dilution of the Osteopontin antibody should be determined by the researcher.

Immunogen

Amino acids DPKSKEEDKHLK were used as the immunogen for this Osteopontin antibody. **Storage** Aliquot and store the Osteopontin antibody at -20oC.