


Osteopontin Antibody (R33573)

Catalog No.	Formulation	Size
R33573-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6696
Localization	Cytoplasmic
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Osteopontin antibody is available for research use only.



Western blot testing of MOLT4 lysate with Osteopontin antibody at 0.5ug/ml. Observed molecular weight: 60~65 kDa (with post-translational modifications), ~35 kDa (unmodified).

Description

Additional name(s) for this target protein: Secreted phosphoprotein 1; SPP1

Application Notes

Optimal dilution of the Osteopontin antibody should be determined by the researcher.

Immunogen

Amino acids DPKSKEEDKHLK were used as the immunogen for this Osteopontin antibody.

Storage

Aliquot and store the Osteopontin antibody at -20oC.