


OTUB1 Antibody (R36379)

Catalog No.	Formulation	Size
R36379-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55611
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This OTUB1 antibody is available for research use only.



Western blot testing of mouse brain lysate with OTUB1 antibody at 0.05ug/ml.
Predicted/observed molecular weight: ~30 kDa.

Description

Additional name(s) for this target protein: OTU domain, ubiquitin aldehyde binding 1

Application Notes

Optimal dilution of the OTUB1 antibody should be determined by the researcher.

Immunogen

Amino acids YKEYAEDDNIYQQK were used as the immunogen for this OTUB1 antibody.

Storage

Aliquot and store the OTUB1 antibody at -20oC.