


# OTUB1 Antibody (R36379)

Catalog No.	Formulation	Size
R36379-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55611
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This OTUB1 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: OTU domain, ubiquitin aldehyde binding 1

## **Application Notes**

Optimal dilution of the OTUB1 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids YKEYAEDDNIYQQK were used as the immunogen for this OTUB1 antibody. **Storage** Aliquot and store the OTUB1 antibody at -20oC.