


OTUB2 Antibody (R36382)

Catalog No.	Formulation	Size
R36382-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	78990
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This OTUB2 antibody is available for research use only.



Western blot testing of human cerebellum lysate with OTUB2 antibody at 1ug/ml.
Predicted/observed molecular weight: ~27kDa.

Description

Additional name(s) for this target protein: OTU domain, ubiquitin aldehyde binding 2

Application Notes

Optimal dilution of the OTUB2 antibody should be determined by the researcher.

Immunogen

Amino acids YKTSHYNILYAADKH were used as the immunogen for this OTUB2 antibody.

Storage

Aliquot and store the OTUB2 antibody at -20oC.