


OTX2 Anitbody [clone 1313CT817.185.162] (F52477)

Catalog No.	Formulation	Size
F52477-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52477-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, k
Clone Name	1313CT817.185.162
Purity	Purified
UniProt	P32243
Applications	Western Blot : 1:2000
Limitations	This OTX2 Anitbody is available for research use only.



Western blot analysis of lysate from Y79 cell line using OTX2 Anitbody at 1:1000.

Description

Probably plays a role in the development of the brain and the sense organs. Can bind to the BCD target sequence (BTS): 5'-TCTAATCCC-3'.

Application Notes

Titration of the OTX2 Anitbody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This OTX2 antibody was produced from a mouse immunized with a recombination protein from human protein.

Storage

Aliquot the OTX2 Antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.