


Ovo-like 2 Antibody / OVOL2 [clone OVOL2/6776] (V4816)

Catalog No.	Formulation	Size
V4816-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V4816-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V4816SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2, kappa
Clone Name	OVOL2/6776
Purity	Protein A/G affinity
UniProt	Q9BRP0
Localization	Nucleus
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT
Limitations	This Ovo-like 2 antibody is available for research use only.



IHC staining of FFPE human colon carcinoma tissue with OVOL2 antibody (clone OVOL2/6776). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.

Description

Ovo-like 2 is a zinc-finger transcription repressor factor.

Application Notes

Optimal dilution of the Ovo-like 2 antibody should be determined by the researcher.

Immunogen

A recombinant partial protein sequence (within amino acids 175-275) from the human protein was used as the immunogen for the Ovo-like 2 antibody.

Storage

Aliquot the Ovo-like 2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.