


# **OXER1 Antibody / 5-oxo-ETE GPCR (RQ7925)**


Catalog No.	Formulation	Size
RQ7925	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

# **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q8TDS5
Applications	Western Blot : 0.5-1ug/ml Flow Cytometry : 1-3ug/million cells Direct ELISA : 0.1-0.5ug/ml
Limitations	This OXER1 antibody is available for research use only.



Flow cytometry testing of human U937 cells with OXER1 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= OXER1 antibody.



Western blot testing of human ThP-1 cell lysate with OXER1 antibody. Predicted molecular weight ~41 kDa.

## **Description**

Oxoeicosanoid receptor 1 (OXER1) also known as G-protein coupled receptor 170 (GPR170) is a protein that in humans is encoded by the OXER1 gene located on human chromosome 2p21; it is the principal receptor for the 5-Hydroxyicosatetraenoic acid family of carboxy fatty acid metabolites derived from arachidonic acid. Enables G protein-coupled receptor activity and icosanoid binding activity. Involved in adenylate cyclase-inhibiting G protein-coupled receptor signaling pathway. Predicted to be located in plasma membrane.

#### **Application Notes**

Optimal dilution of the OXER1 antibody should be determined by the researcher.

#### **Immunogen**

E. coli-derived recombinant human protein (amino acids M1-Q422) was used as the immunogen for the OXER1 antibody.

### **Storage**

After reconstitution, the OXER1 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.