


OXER1 Antibody / Oxoeicosanoid receptor 1 (R36055)

Catalog No.	Formulation	Size
R36055-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	165140
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This OXER1 antibody is available for research use only.



Western blot testing of human lung lysate with OXER1 antibody at 1ug/ml. Predicted molecular weight ~41 kDa.

Description

Additional name(s) for this target protein: Oxoeicosanoid (OXE) receptor 1; 5-oxo-ETE GPCR

Application Notes

Optimal dilution of the OXER1 antibody should be determined by the researcher.

Immunogen

Amino acids SRQWRYREASRK were used as the immunogen for this OXER1 antibody.

Storage

Aliquot and store the OXER1 antibody at -20oC.