


## P130Cas Antibody (R34605)


Catalog No.	Formulation	Size
R34605-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9564
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This P130Cas antibody is available for research use only.



Western blot testing of MCF7 [A] and human breast cancer [B] lysates with P130Cas antibody at 1ug/ml. Predicted molecular weight: ~98kDa.



IHC staining of FFPE human uterus with P130Cas antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## Description

Additional name(s) for this target protein: Breast cancer estrogen resistance 1; BCAR1

## Application Notes

Optimal dilution of the P130Cas antibody should be determined by the researcher.

1. This P130Cas antibody will detect isoforms 1/2/3/5/6/7.

## Immunogen

Amino acids KQTSTFSKQTPHH were used as the immunogen for this P130Cas antibody.

## Storage

Aliquot and store the P130Cas antibody at -20oC.