


p18 Antibody [EEF1E1] (R35937)

Catalog No.	Formulation	Size
R35937-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Pig
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9521
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This p18 antibody is available for research use only.



Western blot testing of K562 lysate with p18 antibody at 0.1ug/ml. Predicted molecular weight: ~20kDa, routinely observed at 18~20kDa.



Western blot testing of pig testis lysate with p18 antibody at 0.1ug/ml. Predicted molecular weight: ~20kDa, routinely observed at 18~20kDa.

Description

Additional name(s) for this target protein: Eukaryotic translation elongation factor 1 epsilon 1 (EEF1E1); ARS-interacting multifunctional protein 3 (AIMP3)

Application Notes

Optimal dilution of the p18 antibody should be determined by the researcher.

Immunogen

Amino acids SKGNKYSAQGERQ were used as the immunogen for this p18 antibody.

Storage

Aliquot and store the p18 antibody at -20oC.