


## p18 Antibody [EEF1E1] (R35937)

Catalog No.	Formulation	Size
R35937-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Pig
<b>Predicted Reactivity</b>	Cow, Dog
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	9521
<b>Applications</b>	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This p18 antibody is available for research use only.



Western blot testing of K562 lysate with p18 antibody at 0.1ug/ml. Predicted molecular weight: ~20kDa, routinely observed at 18~20kDa.



Western blot testing of pig testis lysate with p18 antibody at 0.1ug/ml. Predicted molecular weight: ~20kDa, routinely observed at 18~20kDa.

## Description

Additional name(s) for this target protein: Eukaryotic translation elongation factor 1 epsilon 1 (EEF1E1); ARS-interacting multifunctional protein 3 (AIMP3)

## Application Notes

Optimal dilution of the p18 antibody should be determined by the researcher.

## Immunogen

Amino acids SKGNKYSAQGERQ were used as the immunogen for this p18 antibody.

## Storage

Aliquot and store the p18 antibody at -20oC.