


p27 Antibody (R33948)

Catalog No.	Formulation	Size
R33948-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1027
Applications	Western Blot : 0.5-2ug/ml Immunoprecipitation : suitable under native conditions ELISA (peptide) LOD : 1:32000
Limitations	This p27 antibody is available for research use only.



Western blot testing of human breast lysate with p27 antibody at 0.5ug/ml. Predicted molecular weight ~27kDa.

Description

Additional name(s) for this target protein: P27KIP1, CDKN1B

Application Notes

Optimal dilution of the p27 antibody should be determined by the researcher.

Immunogen

Amino acids EQTPKKPGLRRRQT were used as the immunogen for this p27 antibody.

Storage

Aliquot and store the p27 antibody at -20oC.