


## p27 Antibody (R33948)

Catalog No.	Formulation	Size
R33948-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Cow, Dog, Mouse, Rat
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	1027
<b>Applications</b>	Western Blot : 0.5-2ug/ml Immunoprecipitation : suitable under native conditions ELISA (peptide) LOD : 1:32000
<b>Limitations</b>	This p27 antibody is available for research use only.



Western blot testing of human breast lysate with p27 antibody at 0.5ug/ml. Predicted molecular weight ~27kDa.

## Description

Additional name(s) for this target protein: P27KIP1, CDKN1B

## Application Notes

Optimal dilution of the p27 antibody should be determined by the researcher.

## **Immunogen**

Amino acids EQTPKKPGLRRRQT were used as the immunogen for this p27 antibody.

## **Storage**

Aliquot and store the p27 antibody at -20oC.