


p53 Antibody / TP53 [clone RM387] (R20404)

Catalog No.	Formulation	Size
R20404-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM387
Purity	Protein A purified from animal origin-free supernatant
UniProt	P04637
Localization	Nuclear
Applications	Immunohistochemistry (FFPE) : 1:100-1:200 Western Blot : 1:1000-1:2000
Limitations	This recombinant p53 antibody is available for research use only.



Western blot testing of human HEK293 cell lysate with recombinant p53 antibody at 1:1000 dilution. Expected molecular weight ~53 kDa.



IHC staining of FFPE human lung cancer tissue with recombinant p53 antibody at 1:100.

Description

This antibody reacts to human p53 / Tumor protein 53 / TP53.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant p53 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to human p53 was used as the immunogen for the recombinant p53 antibody.

Storage

Store the recombinant p53 antibody at -20oC.