


P70 S6 Kinase Antibody / P70S6K [clone RM438] (R20453)

Catalog No.	Formulation	Size
R20453-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM438
Purity	Protein A purified from animal origin-free supernatant
UniProt	P23443
Applications	Immunohistochemistry (FFPE) : 1:100-1:400 Western Blot : 1:1000-1:2000
Limitations	This recombinant P70 S6 Kinase antibody is available for research use only.



IHC staining of FFPE human brain cancer tissue with recombinant P70 S6 Kinase antibody at 1:100.



IHC staining of FFPE human breast cancer tissue with recombinant P70 S6 Kinase antibody at 1:400.



Western blot testing of human MCF7 cell lysate with recombinant P70 S6 Kinase antibody at 1:2000 dilution. Predicted molecular weight: 60-70 kDa.

Description

This antibody reacts to P70 S6 Kinase/p70S6K.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant P70 S6 Kinase antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the N-terminus of p70S6K was used as the immunogen for the recombinant P70 S6 Kinase antibody.

Storage

Store the recombinant P70 S6 Kinase antibody at -20°C.